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Instructions for use
Pre-desinfection - Cleaning - Sterilization
of reusable medical devices

Reusable surgical instruments - Reusable sterilization baskets and accessories

General recommendations	Do not use this product for any purpose other than that for which it is intended. This device is intended for use by trained and qualified healthcare professionals. It may not be used for purposes other than those for which it is intended. Any serious incident in connection with the device should be notified to the manufacturer and the competent authority of the member State in which the user and/or patient is established.
Recommendations for choosing detergent solutions	For STAINLESS STEEL devices : choose alkaline cleaning solutions, sodium hypochloride (bleach) is totally contraindicated (causes corrosion). Avoid ALL CONTACT with corrosive substances such as chlorine, chlorides, iodine and fluoride. For ALUMINIUM ALLOY devices: - Use detergent with a neutral pH (pH = 7-8) - Do not use sodium hydroxide-based products - Pion inactivating detergent may be highly alkaline (pH> 8) and therefore not suitable to clean anodized aluminum devices without the use of a suitable neutralizing agent; disinfection and final rinse shall be performed with, at least, a demineralised water. Do not use products containing chlorine or chloride. Follow the instructions for use and the manufacturer's recommendations closely.
Recommendations for water quality	The quality of the water used for the various device retreatment operations is very important: it must be evaluated, controlled and monitored. Inappropriate water chemical composition will have a negative impact on both the treatment process and the external appearance of the devices. Avoid hard water and ensure that chlorine, calcium carbonate and iron contents are low. Demineralised water is recommended for the final rinsing step in all cases.
Limits of reconditioning	Repeated sterilization cycles have little effect on the devices. The lifespan of a device is dictated by wear and damage related to use. It is important to ensure that the various detergents used during reconditioning are compatible.
Before first use of new devices	Visually inspect the surgical instruments to ensure they are in good working order. Clean and disinfect several times to ensure that any eventual residues are removed before they are first sterilised and used.

Instructions

Pre-disinfection	The devices should be reconditioned as soon as possible after use. The devices should be pre-disinfected to reduce the micro-organism count and facilitate subsequent cleaning. Organic material must not be allowed to dry on the devices. Pre-disinfection must be done as soon as possible after the medical device is used and as close as possible to the place of use, prior to cleaning. Use a neutral pH pre-disinfection solution following the manufacturer's instructions for quantities/concentrations and soaking times. It is preferable to use water with a HT between 4° et 8° F. The oxidant chlorine content must not be more than 100 mg/l. The risk of corrosion is high if the content exceed 100 mg/l. Hinged surgical instruments must be widely opened in order to access to every areas. All surfaces, cavities, lumens and openings must be covered with the pre-disinfection solution. After pre-disinfection, rinse the device thoroughly with preferably soft or demineralized water.
Preparation for cleaning	Remove excess organic material with a disposable cloth / paper. Open the instruments. Surgical instruments formed of several parts must be dismantled. The inner cannula must be removed from cannulas with a small diameter for cleaning. Identify and discard any oxidised device.
Automated cleaning in a washing disinfecting machine	Prepare and load the instruments into the machine in compliance with the type of equipment used. The instruments must be positioned so that all joints will be washed but not damaged. Select the appropriate program for the type of load in compliance with the manufacturer's instructions of the automatic cleaning machine. It is better to use a detergent with a neutral pH (7-8). Do not use an acid detergent (pH<7). Detergents containing chlorine are not appropriate. Follow the instructions for concentration, temperature and contact time given by the manufacturer of the cleansing agent. Use demineralised water or water purified by reverse osmosis for the final rinsing step.
Manual cleaning	Manual cleaning involves brushing and/or soaking the instruments in a detergent and disinfectant solution. After cleaning, the instruments must be rinsed carefully and thoroughly with tap water. Use demineralised water or water purified by reverse osmosis for the final rinsing step. Do not use metallic brushes or sponges or abrasive cleaning products. It is preferable to use a detergent with a neutral pH. Do not use an acid detergent (pH<7). Follow the instructions for concentration, mechanical action, temperature and contact time provided by the manufacturer of the cleansing agent. Make sure that hinged surgical instruments are cleaned in both open and closed positions. Clean cannulas and holes with an appropriate brush, ensuring that the brush reaches the bottom of the instrument. During rinsing, make sure that the water runs through the cannulas and that blind holes are filled and emptied several times.

Drying	Automated cleaning in a washing-disinfecting machine: - the drying phase must be sufficiently long,- do not exceed a temperature of 120°C Manual drying: - dry the instrument with a clean, absorbent, soft, lint-free cloth, - dry cannulas and lumens with a medical air system
Maintenance	Mobile parts (hinges and closures) must be regularly lubricated with a medical, sterilisable, steam-permeable maintenance oil. Check the general and surface conditions of the device and discard any device blunt, damaged, defective or with deep marks and/or scratches..

Controls and tests	Allow the instrument to cool at room temperature. Lubricate mobile parts with a medical, sterilisable, steam-permeable maintenance oil. After each instrument cleaning operation, check the following: - all traces of organic material have been removed: if the instrument is not sufficiently clean, repeat the cleaning phase; - the instrument is in good working order and has not been damaged (no bent, broken, cracked, worn or split parts and cannula/trocar openings not blocked); cutting edges must be continuous and free of nicks; - hinged instruments: check that the hinges move easily but are not loose;- verify the locking mechanisms (rack systems, hooks); - assemble any instruments that can be dismantled and check that they are working correctly; when the instruments are made of several pieces, fit them together with the other parts. - Check cutting edge blades. Refer to the table below for scissors:																		
	<table><tr><th>Materials</th><th>Number of thicknesses</th><th>Instrument category</th><th>Examples</th></tr><tr><td>Latex (e.g. surgical gloves)</td><td>1</td><td>Microsurgery spring scissors</td><td>Barraquer scissors Iridectomy scissors Wecker scissors</td></tr><tr><td rowspan="3">Gauze No. 18</td><td>2</td><td>Light with rings</td><td>Iris scissors Strabismus scissors Stevens tendon scissors Kilner dissecting scissors Potts de Martell , Joseph scissors</td></tr><tr><td>5</td><td>Middle with rings</td><td>Mayo, Metzenbaum scissors Nail scissors Sims uterine scissors Dauphin bandage scissors Tonsil scissors</td></tr><tr><td>8</td><td>Heavy</td><td>Lister bandage type scissors Smith bandage scissors Intestine scissors Doyen tissue scissors</td></tr></table>	Materials	Number of thicknesses	Instrument category	Examples	Latex (e.g. surgical gloves)	1	Microsurgery spring scissors	Barraquer scissors Iridectomy scissors Wecker scissors	Gauze No. 18	2	Light with rings	Iris scissors Strabismus scissors Stevens tendon scissors Kilner dissecting scissors Potts de Martell , Joseph scissors	5	Middle with rings	Mayo, Metzenbaum scissors Nail scissors Sims uterine scissors Dauphin bandage scissors Tonsil scissors	8	Heavy	Lister bandage type scissors Smith bandage scissors Intestine scissors Doyen tissue scissors
	Materials	Number of thicknesses	Instrument category	Examples															
	Latex (e.g. surgical gloves)	1	Microsurgery spring scissors	Barraquer scissors Iridectomy scissors Wecker scissors															
	Gauze No. 18	2	Light with rings	Iris scissors Strabismus scissors Stevens tendon scissors Kilner dissecting scissors Potts de Martell , Joseph scissors															
5		Middle with rings	Mayo, Metzenbaum scissors Nail scissors Sims uterine scissors Dauphin bandage scissors Tonsil scissors																
8		Heavy	Lister bandage type scissors Smith bandage scissors Intestine scissors Doyen tissue scissors																
Following the instructions in the table, the materials used for the test must be cut in one go and three times using two-thirds of the blade without exerting lateral pressure. It must be cut smoothly, cleanly and without tearing.																			
Packaging	Protect the ends of fragile or sharp instruments with suitable perforated silicone protectors that are able to withstand the sterilization conditions and are permeable to the sterilizing agent. Place device with fine ends and/or microsurgery instruments on appropriate stands (e.g. silicone rests). Hinged instruments should be set to the first notch. Packaging is performed in compliance with the ISO 11607-1 & 2 standards.																		
Sterilization	Sterilization using a validated steam sterilization process. Autoclave at 134°C for at least 18 minutes. The water used to produce the steam must be compatible with the sterilization process. The sterilization process is subject to national legislation and guidelines and it is therefore impossible to describe a single, universally-applicable procedure. It is the user's responsibility to ensure that the expected results are actually obtained with the sterilization process used (methods, equipment, materials and staff). The user is therefore responsible for validation and routine monitoring of all the processes concerned.																		
Storage	The device must be packaged, identified and stored protected from direct sunlight, humidity and all sources of contamination.																		
Further information	The device can be severely damaged by acids in the products used to recondition them. It is preferable to return devices to MEDLANE for reconditioning. Written proof must be provided showing that instruments returned for maintenance and/or repair have undergone all the procedures applicable for sterile medical devices, both prior to delivery and before they are sent back. MEDLANE is not responsible for any incidents or accidents caused by devices that are dismantled, transformed or repaired by a third party unauthorised by MEDLANE. Devices whose lifespan has expired and / or are no longer able to function as the intended use must be discarded properly according to internal procedures defined by the health facility.																		

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The instructions supplied above have been approved by the medical device manufacturer for the preparation of medical devices for reuse. It is the competent department's responsibility to ensure that the Sterilization process actually used and involving the equipment, materials and staff required for reconditioning, achieves the expected result. This generally requires validation and routine controls on the processes used. It is also important to analyse any deviations from the manufacturer's instructions for Sterilization to establish whether the process is still effective or whether such deviations have had a negative impact.